INSTRUCTIONS FOR USE

Intramedullary Nailing Implants

Document 700001 | Revision 17

Select a language to open the corresponding Instructions for Use.

[English](#)

EN

[Czech](#)

CS

[Italian](#)

IT

[Romanian](#)

RO

INSTRUCTIONS FOR USE / WARNINGS & PRECAUTIONS

AUSTOFIX INTRAMEDULLARY NAILING SYSTEM

Intramedullary Nailing Implants

Associated device systems with these instructions for use

- Austofix F1 Proximal Femoral Nail
- Austofix F3 Proximal Femoral Nail
- Austofix S2 Supracondylar Nail
- Austofix UTN Tibial Nail
- Austofix CHN Humeral Nail


GENERAL DESCRIPTION

The Austofix Intramedullary Nailing System is a family of titanium implants used for fracture fixation in orthopaedic trauma surgery. The Austofix nailing system is implanted surgically by trained orthopaedic surgeons to fix fractures in any of the three long bones: humerus, femur and tibia. There are two distinct types of implants in this system. Intramedullary nails are long narrow solid or tubular rods of circular cross-section, with cross holes for placement of locking screws. They are designed to be implanted inside the hollow 'medullary canal' of the humerus, femur and tibia. Intramedullary screws are bone screws, which fix the nail in position relative to the bone and prevent rotation. The nail and screws function as a type of internal splint, holding the broken bone in place until it grows together again. Around the one-two month period, the fracture should have healed sufficiently to be taking the full weight of the patient, and the implants no longer perform any function. It is safe for them to remain in situ indefinitely, but they are often removed approximately one year after insertion. Hence, the implants themselves have an indefinite life either before or after implantation.

The implants are available in different configurations and sizes for different fracture types and anatomies. However, the same basic fixation techniques apply for each type. The instrument set includes a basic set of surgical instruments, K-wires, drills, and specific drill guide instruments to ensure accurate alignment of the implants. Each instrument set is clearly marked with the product name.

All implantable devices are single use only. Nails and screws are supplied sterile (gamma irradiation) in two layers of OPA packaging. All intramedullary nails and screws are manufactured from implant grade titanium alloy in accordance with ISO 5832-3.

DEFINITION OF SYMBOLS

	Manufacturer		Sterilised using irradiation		Consult IFU or electronic IFU
	Authorised representative in the European Union		MR Conditional		Caution
	Date of manufacture		Do not use if package is damaged, consult IFU		Medical device
	Use-by date		Do not re-use		Unique device identifier
	Batch code		Do not resterilise		
	Catalogue number		Double <i>sterile</i> barrier system		

INTENDED PURPOSE

Intramedullary Nailing Implants are intended to be used for temporary fixation and stabilization of long bones in various anatomical regions such as proximal femur, femoral shaft, tibia and humerus.

INDICATIONS AND CONTRAINDICATIONS

Indications

Indications for interlocking intramedullary nails include severely comminuted, spiral, long oblique and segmental fractures, non-unions and malunions, and bone lengthening/shortening. The general principles of patient selection and sound surgical judgement apply. The size and shape of the long bones present limiting restrictions on the size and strength of implants. Specific indications for intramedullary nailing depend on the anatomical location. Devices are designed to be used only by trained orthopaedic surgeons in a hospital environment.

Contraindications

The general principles of patient selection and sound surgical judgement apply. Allergies and other reactions to device materials, although infrequent, should be considered, tested for (if appropriate), and ruled out preoperatively. Contraindications to be avoided include:

- Patients with open epiphyseal plates.
- Insufficient quantity or quality of bone, conditions which tend to retard healing, and blood supply limitations.
- Previous or active infection.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Conditions that tend to affect the patient's ability or willingness to restrict activities during the healing period.

Skeletal deformity precluding nail use or obliterated medullary canal.

INTENDED USER

Intramedullary Nailing Implants are intended for use by qualified orthopaedic surgeons trained in intramedullary nailing and fluoroscopic imaging for accurate placement and fixation.

INSTRUCTIONS FOR USE / WARNINGS & PRECAUTIONS

AUSTOFIX INTRAMEDULLARY NAILING SYSTEM

PATIENT TARGET GROUP

Adult patients with long bone fractures of the femur, tibia, or humerus. These devices are particularly suitable for individuals with good bone quality but may also be used in patients with osteoporosis, provided adequate fracture stability can be achieved.

WARNINGS AND PRECAUTIONS

Preoperative

1. Use care in handling and storage of implant components. Cuffing, sharply bending, or scratching the surface can significantly reduce the strength and fatigue resistance of the implant system. This, in turn, could induce cracks and/or invisible internal stresses that could lead to fracture of the implants. Implants and in storage should be protected from corrosive environments such as salt air, moisture, etc.
2. Patient conditions and/or predispositions, such as those addressed in Contraindications above, should be avoided.
3. An adequate inventory of implant sizes should be available at the time of surgery.
4. Allergies and other reactions to device materials, although infrequent, should be considered, tested for (if appropriate), and ruled out pre-operatively.
5. Before the initial use of these implants, we recommend that the surgeon acquaint himself with them and attend a technique seminar. Surgical Technique brochures are available upon request at no charge and should be reviewed by the surgeon prior to initial surgery. Skill in the use of this technique should be acquired on less complicated fractures before attempting its use in unstable, difficult fractures. As a general guide, reaming to a diameter at least 1.5mm greater than the nail should always be considered.
6. The patient should be advised that a second more minor procedure for the removal of implants may be necessary.

Operative

1. Selection of the proper nail length and diameter is extremely important and must be carefully sized to the patient, taking into account the patient's age, weight, and cortical bone quantity. As a rule, the largest implant that easily fits the canal should be used. Small canals require enlargement by reaming.
2. Inspection and trial assembly are recommended prior to implantation to determine if instrument components or implants have been damaged during storage or prior procedure.
3. Care should be taken not to scratch, bend sharply, or cut metal components during surgery for the reasons stated.
4. Refer to the outer carton labels, surgical technique, or product catalogue for information on the correct size of screws for each nail.
5. A stable construct should be achieved and verified by Xray imaging.
6. Once removed, implants should never be reused since internal stresses (in the implant) that are not visible may lead to early bending or fracture.
7. Excessive drilling or reuse of drills can produce drill wear, bluntness and heat generation, leading to increased operating time and potential osteonecrosis.

Implantable devices labelled with "Do not resterilise" symbol must not be re-sterilised. Re-sterilisation of implantable devices can result in product not being sterile and/or not meeting performance specifications and/or altered material properties.

Postoperative

1. Although Austofix Nails are designed for maximum strength and performance, it must be well understood that intramedullary nails are not intended to carry the full load of the patient acutely nor for extended periods of time. All patients should be cautioned against significant weight bearing prior to good callus formation. For this reason, patients who are obese and/or non-compliant, as well as patients who could be predisposed to delayed or non-union, must have auxiliary support. The implant may be exchanged for a larger, stronger nail subsequent to the management of soft tissue injuries.
2. Additional postoperative precautions should be taken when the fracture line occurs within 5cm of the nail's screw hole, as this situation places greater stress on the nail.
3. Postoperative directions and warnings to patients by physicians, and appropriate nursing care, are extremely important, particularly those admonitions that concern early weight-bearing or active use of the extremities. These activities substantially increase the stress on implants that can lead to complications.
4. Periodic X-ray examinations for at least the first three (3) months postoperatively are necessary to detect changes in position, non-union, loosening, bending, or cracking of components. With evidence of these conditions, patients should be closely observed, the possibilities of further deterioration evaluated, and the benefits of reduced activity and early revision considered.
5. Early weight bearing should be considered only in those cases with stable fractures and good bone-to-bone contact.

INFORMING THE PATIENT

The surgeon is responsible for informing the patient of the expected benefits, residual risks, contraindications, and possible undesirable side-effects associated with the Austofix Intramedullary Nailing System. The patient should be informed that the implant does not replicate normal bone, may fail if subjected to excessive load or non-compliance, and that additional procedures, including implant removal, may be required.

MRI Safety Information

Austofix devices should be considered MR Conditional. Devices have minimal ferro-magnetism with minimal risk in strong magnetic fields, since devices are fixed in bone. Austofix recommends the following when using INS in an MRI environment.

- Static magnetic field strength (T) 1.5T or 3.0T; Maximum special field gradient 30T/m (3000 gauss/cm); RF excitation Circularity Polarised (CP); RF transmit coil type Integrated Whole Body Transmit Coil; Operating mode Normal Operating Mode; Maximum whole-body SAR (W/kg) 2 W/kg (Normal Operating Mode).
- Scan Duration – 1.5T 2 W/kg whole-body average SAR for 6 minutes of continuous RF (a sequence or back-to-back series/scan without breaks) followed by a wait time of 6 minutes if this limit is reached, for the total scanning session duration of up to 1 hour (or 60 minutes); 3.0T 2 W/kg whole-body average SAR for 1 hour (or 60 minutes) of continuous RF (a sequence or back to back series/scan without breaks).
- MRI image artifact - The presence of this implant can produce an image artifact of approximately 27 mm from the system when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

CLINICAL BENEFITS

The Austofix intramedullary nailing system provides stable internal fixation for long bone fractures, enabling early mobilization and reduced complication rates compared to conventional fixation methods.

PERFORMANCE CHARACTERISTICS

When used according to the manufacturer's Instructions for Use and Labelling, the intermedullary nailing implants are intended to function as a type of internal splint, holding the broken bone in place to support healing until bone fusion has been achieved.

EXPECTED LIFETIME

The nail and screws function as a type of internal splint, holding the broken bone in place until it grows together again. Around the one-to-two-month period, the fracture should have healed sufficiently to be taking the full weight of the patient, and the implants no longer perform any function. It is safe for the device to remain in situ indefinitely. This depends on the patient situation and indication, but they are often removed approximately one year after insertion.

INSTRUCTIONS FOR USE / WARNINGS & PRECAUTIONS

AUSTOFIX INTRAMEDULLARY NAILING SYSTEM

POSSIBLE COMPLICATIONS AND SIDE EFFECTS

1. Loosening, bending, cracking, or fracture of the nails or screws, or loss of fixation in the bone, attributable to the factors listed in Contraindications above and/or Warnings and Precautions below.
2. Loss of anatomic position with non-union or malunion with rotation or angulation.
3. Infections, both deep and superficial.
4. Fat embolism syndrome.
5. Allergies and other reactions to device materials.
6. Irritation of soft tissues, including impingement syndrome.

In the event of a serious incident involving an Austofix product, users must contact the manufacturer and the Competent Authority of the relevant Member State.

PACKAGING AND LABELLING

Implants labelled as sterile have been sterilised by a minimum of 25 kGy of gamma irradiation. Inspect packaging for punctures or other damage prior to surgery. All implants that are provided sterile should be accepted only if the factory packaging and labelling arrive intact. If the sterile barrier has been compromised in any way, the devices should not be used. Any such instances should be reported to the manufacturer and the devices returned via the supplier for evaluation by the manufacturer.

RESTERILISATION

Only the reusable instrument sets are intended to be sterilised by the user, and may be resterilised, if necessary, by steam autoclaving in appropriate protective wrapping. The following process parameters are recommended for these devices: Manual Clean: 10L water (max temperature 35°C) dosed with 4mL/L of Viruzyme V. Then Automatic clean: Prewash cold water (3 min), Wash: Neutral Enzymatic Chemistry (Amity Viruzyme-V) dosed at 4mL/L (192ml) @ 60°C, 12 min, Rinse: Incoming Hot Water, 1 min, Thermal Disinfection: Reverse Osmosis Water, 93°C @ 1 min, 0.3mL/L (14, 4ml) Rinse Aid pH Neutral. Sterilise: Pre-vacuum cycle, 4 minutes at 134 C, followed by 60 minutes of drying time. Detailed instructions for reprocessing, including cleaning, disinfection and resterilisation, are provided in the companion document: *F40-LG-07 General Requirements for Reprocessing*.

DISPOSAL

Hospitals should dispose of explanted devices using established procedures for clinical waste.

The Summary of Safety and Clinical Performance (SSCP) will be available on the European Database on Medical Devices (EUDAMED) once the module becomes functional.

Australian Orthopaedic Fixations Pty Ltd
(t/a Austofix)
28 Dalglish Street, Thebarton, SA 5031
AUSTRALIA
Ph: +61 8 8351 0644
info@austofix.com.au
www.austofix.com.au

Advena Ltd. Tower Business Centre, 2nd Flr.,
Tower Street, Swatar,
BKR 4013 Malta

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