




Flexible Reamer

PRODUCT BROCHURE

Austofix Flexible Reamer


The Austofix Flexible Reamer Set includes a 470 mm flexible reamer shaft manufactured from hardened stainless steel, offering durability and resistance to wear. Reamer heads are available in 0.5 mm increments ranging from Ø8 mm to Ø14 mm, providing precise canal preparation.

An optional 210 mm extension is available for suprapatellar access, specifically designed to interface with the Austofix Universal Tibial Nail.



Quick-Connect Locking Mechanism

- Simple push-and-twist action for rapid assembly and removal
- Secure locking interface ensures stable transmission of torque during reaming
- Minimises surgical delay with fast head transfer and reduces handling errors with reliable engagement




Flexible Reamer Shaft 543047

- 470mm Shaft for femoral, humeral and tibial nail insertion
- Z/H Connection
- Sterile Consumable Item


Reamer Extension 543046

- 210mm Extension specific for tibial suprapatellar approach
- Quick-Connect Locking with 470mm Shaft
- Sterile Consumable Item




Removal Caddy 621052

- Designed for the safe and efficient removal of used reamer heads
- Quick-release mechanism allows simple slot-and-pull detachment
- Protects surgical gloves from damage




Guidewire Gripper 543049

- Secure grip to advance the guidewire
- Compatible with multiple wire diameters
- Durable Stainless-Steel construction



T-Handle 531004


- Manual, reliable control for hand reaming
- Z/H Coupling with Silicon Handle
- Ø4mm Cannulation



- | | | | | | |
|----|--------|---------------------------------|-----|--------|---------------|
| 1. | 543047 | Flexible Reamer Shaft - Sterile | 10. | 543105 | Ø10.5mm Head |
| 2. | 543046 | Extension Shaft - Sterile | 11. | 543110 | Ø11mm Head |
| 3. | 543049 | Guidewire Gripper | 12. | 543115 | Ø11.5mm Head |
| 4. | 531004 | T-Handle | 13. | 543120 | Ø12mm Head |
| 5. | 543008 | Ø8mm Head - End Cutting | 14. | 543125 | Ø12.5mm Head |
| 6. | 543085 | Ø8.5mm Head | 15. | 543130 | Ø13mm Head |
| 7. | 543090 | Ø9mm Head | 16. | 543135 | Ø13.5mm Head |
| 8. | 543095 | Ø9.5mm Head | 17. | 543140 | Ø14mm Head |
| 9. | 543100 | Ø10mm Head | 18. | 621052 | Removal Caddy |



Austofix



Legal Manufacturer

Australian Orthopaedic Fixations Pty Ltd
18 Kinkaid Avenue, North Plympton SA 5037, AUSTRALIA

Telephone: 1300 727 075
Fax: 1300 727 380
Email: sales@austofix.com.au



Advena Ltd, Tower Business
Centre, 2nd Flr, Tower Street,
Swatara, BKR 4013 Malta

DISTRIBUTOR STICKER

700409 Flexible Reamer Brochure Rev 0

Australian Made • Australian Owned • Australian Designed

austofix.com.au